**Inhalt:**

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**Contents:**

- |  |                                 |
|--|---------------------------------|
| 1 Aluminium container with 100 test sticks |                                 |
| 1 Bottle Aluminium-1                       | 1 Measuring tube with 5 ml mark |
| 1 Bottle Aluminium-2                       | 1 Test tube                     |

**Safety precautions:**

**Aluminium-1 contains 50% potassium hydroxide solution. Causes severe burns. Aluminium-2 contains 24% acetic acid. Irritating to eyes and skin.**

Keep locked up and out of reach of children. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable protective gloves and eye protection. If you feel unwell, seek medical advice immediately.

**Method of application:**

1. Rinse measuring tube with the sample and fill to the mark.
2. Add six drops of Aluminium-1 (potassium hydroxide solution) and shake carefully (pH value 13.5 - 14).
3. Remove only as many test sticks as are required, and reseal the container immediately after use. Do not touch the test paper zone.
4. Dip the test stick briefly into the prepared test solution and shake of excess liquid.
5. Place test tube into the cavity of the thermoformed lining and fill with six drops Aluminium-2.
6. Place test stick into the filled test tube. After 1-2 minutes remove test stick and compare the test paper zone with the colour scale. In the presence of aluminium the test paper turns red (pink = 0 mg/l).

A red colouration of the unused test paper is normal and does **not** indicate decomposition of the reagents.

**Interferences:**

Beryllium(II) interferes. Copper(II) interferes in excess of 10 mg/l (can be precipitated with KI or cadmium powder). The following anions cause low test results:

- > 100 mg/l  $MnO_4^-$ ,  $PO_4^{3-}$ ,  $S^{2-}$ ,  $SO_3^{2-}$ .
- > 500 mg/l  $B_4O_7^{2-}$ ,  $F^-$ ,  $Fe(CN)_6^{3-}$ ,  $S_2O_4^{2-}$ ,  $S_2O_5^{2-}$ ,  $SiO_4^{4-}$ , citrate and tartrate.

**Storage:**

Avoid exposing the sticks to sunlight and moisture. Store the kit below +30 °C in a dry place.

MACHEREY-NAGEL GmbH & Co. KG · P.O. Box 10 13 52 · D-52313 Dueren (Germany)  
 Phone +49 24 21 9 69-0 · Fax +49 24 21 9 69-199 · E-Mail: sales@macherey-nagel.de

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