QUANTOFIX® Cyanide



Contents:

1 Aluminium container with 100 test sticks

1 Bottle Cyanide-1 1 Measuring tube with 5 ml mark

1 Bottle Cyanide-2 1 Small measuring spoon

Safety precautions:

Keep out of reach of children. When performing the test, do not eat, drink or smoke. When you have finished the test, flush the sample with plenty of water in a place, where contact with food or dishes is not possible.

Method of application:

- 1. Rinse measuring tube with the sample and fill to the 5 ml mark.
- 2. Measure pH value (e.g. with pH-Fix 0-14 Cat. No. 921 10).
- 3. If the pH is >10, add dropwise diluted sulphuric acid, until the pH is 6-7.
- 4. Add one measuring spoon Cyanide-1 (buffer mixture) and shake carefully.
- Remove only as many test sticks as are required, and reseal the container immediately after use.
- 6. Add 5 drops Cyanide-2 to the sample solution and shake carefully.
- 7. **Immediately** place test stick into the sample solution.
- After 45 seconds remove test stick and within 10 seconds compare with the colour scale.

Interferences:

Only free cyanides, and cyano complexes which can be decomposed with chlorine, are covered by this test. Thiocyanates (SCN⁻) above 1 mg/l show a similar coloration. lodide and bromide above 5 mg/l cause a decrease in colour intensity or – in higher concentrations – even inhibit the colour formation. Sulphide ions (S²⁻) up to 20 mg/l do not interfere. Higher concentrations also decrease or inhibit the colour formation.

If the pH is >10, then the capacity of the buffer (Cyanide-1) is not sufficient to adjust the sample to a pH of 6-7. In this case the pH value should be adjusted to a pH of 6-7 by dropwise addition of diluted sulphuric acid prior to the test.

Storage:

Avoid exposing the sticks to sunlight and moisture. Store the kit below +30 $^{\circ}\text{C}$ in a dry place.