Specification

Product:	Chromafil [®] Polytetrafluoroethylene (PTFE)
Cat. No.:	729005
Filters are HPLC certified for low level UV-absorbing substances for gaining	
accurate chromatographic results	
Housing material:	Polypropylene (PP)
Membran filter material:	Polytetrafluoroethylene
Pore diameter:	0.45 μm
Filter diameter:	25 mm
Filtration area:	$\sim 4.9 \text{ cm}^2$
Connections inlet:	Female luer
Connections outlet:	Male luer
Dead volume:	~120 µl

HPLC-test

Conditions: Two ml of the solvent (specified under the chromatograms) were applied to the filter; $100 \,\mu$ l of the filtrate were injected to the HPLC. Eluent A: water, eluent B: acetonitrile Gradient 10% -> 95% B in 25 minutes Flow rate: 0.5 ml/min Sensitivity: -5 to 10 mAU at 254 nm Column: 125/4 NUCLEODUR C18 Gravity 5µm (Cat. No. 760100.40)

www.mn-net.com




Germany and international: Tel.: +49 (0) 24 21 96 90 Fax: +49 (0) 24 21 96 91 99 e-mail: sales-de@mn-net.com


MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6-8 · D-52355 Düren · Germany Switzerland: MACHEREY-NAGEL AG Tel.: +41 (0) 62 388 55 00 Fax: +41 (0) 62 388 55 05 e-mail: sales-ch@mn-net.com


France: MACHEREY-NAGEL EURL Tel.: +33 (0) 3 88 68 22 68 Fax: +33 (0) 3 88 51 76 88 e-mail: sales-fr@mn-net.com

USA: MACHEREY-NAGEL Inc. Tel.: +1 610 559 98 48 Fax: +1 610 559 98 78 e-mail: sales-us@mn-net.com


Specification







Acetonitrile-Water (1:1) filtrate



Switzerland: MACHEREY-NAGEL AG Tel.: +41 (0) 62 388 55 00 Fax: +41 (0) 62 388 55 05

e-mail: sales-ch@mn-net.com

Acetonitrile filtrate

MN

(MN)

Germany

and international: Tel.: +49 (0) 24 21 96 90 Fax: +49 (0) 24 21 96 91 99 e-mail: sales-de@mn-net.com

www.mn-net.com

